

SKYN[®] Erbium Laser Skin Resurfacing

Treatment:	SKYN [®] Erbium Laser Skin Resurfacing or skin resurfacing; laser resurfacing; Action II erbium YAG laser; erbium YAG laser; Dermatology laser; Er:YAG.
Technology:	 <p>Lutronic Action II</p>
What:	Laser resurfacing is typically recommended for clients who wish to reduce fine lines and wrinkles, pigmentation from sun exposure or acne scarring. An erbium YAG laser is used as a tool to resurface the skin to reveal desirable, improved changes in the skin.
Treats:	Lines and wrinkles, collagen stimulation, pigmentation, acne & scarring, skin tightening, sun damage, skin rejuvenation, skin tightening.

Erbium Laser Skin Resurfacing

Laser resurfacing uses ablative lasers to reduce the appearance of wrinkles and scars, even pigmentation, tighten skin and remove lesions.

Treatment with Erbium laser is currently the method "minimally invasive" more effective and tolerated to obtain a resurfacing of the skin. This treatment makes it possible to obtain a net and marked improvement of the treated skin (face, neck, cleavage, hands, recent scars and stretch marks) without excessive discomfort by the patient and short recovery times.

SKYN's erbium laser skin resurfacing with Action II is the treatment of choice for clients who want to achieve dramatic improvement in the appearance of their skin. This treatment successfully removes superficial lines to moderately deep wrinkles with less pain, fewer side-effects, and a more rapid recovery

than other treatments. An erbium laser provides a milder and less invasive treatment option than a carbon dioxide (CO2) laser.

How does Erbium Laser Resurfacing Work?

Erbium laser skin resurfacing precisely ablates the epidermis (the outermost layer of the skin) to reveal the underlying skin that is smoother and younger looking. The procedure also activates collagen to promote further tightening of the skin. Our Doctors will precisely remove the necessary layers of skin to eliminate or reduce superficial to moderately deep lines and wrinkles.

Erbium Laser Skin Resurfacing Candidates

A candidate for erbium laser skin resurfacing is someone who is seeking to correct fine lines and wrinkles or achieve a younger and more refreshed look.

Benefits of Erbium Laser Treatment

Erbium laser skin resurfacing has been found safe and effective for the effective treatment of fine lines and moderate wrinkles. Clients who have previously had laser skin resurfacing to enhance their look and improve their self-confidence now find themselves with greatly improved, smoother, and younger-looking skin. Erbium laser skin resurfacing has also been shown to stimulate the production of collagen, which further smoothens and tightens the skin.

Is Erbium Laser Skin Treatment Safe?

Even as erbium laser skin resurfacing provides a highly effective treatment for fine lines and wrinkles with a shorter recovery period and fewer side effects, those receiving treatment should still expect their skin to appear pink or red after the procedure. Some swelling may also occur.

Before + After



Acne Scar Treatment



Seborrheic Keratosis Treatment

The versatile fractional rejuvenation Er:YAG laser solution

Action II's fractional technology leads to a rapid healing process with minimum downtime. It allows physicians to provide a virtually painless treatment.

Applications

- Skin Rejuvenation
- Acne Scars
- Wrinkles
- Large Pores
- Skin Toning
- Texture Improvement
- Lesions

Testimonials

“Er:YAG laser is a safer device that can remove various epidermal and dermal lesions with minimal damage to surrounding normal tissue. ACTION also has fractional mode and long pulse mode, so we are able to utilize this machine for many cosmetic problems such as photo-damaged skin and warts without bleeding.”

Joo Yeon Ko, MD, South Korea

“ACTION II is a multi-functional 2940 nm Er:YAG laser providing conventional resurfacing, but with the addition of three other modes: the long pulsed, fractional and unique ‘shining peel’ modes ACTION II is a truly outstanding choice for the single-doctor clinic, because the system offers unparalleled ease of use with safety, solidity and breakdown-free reliability.

I can confidently recommend ACTION II as an essential device for the state-of-the-art clinic, representing the most ultra-advanced Er:YAG laser in the world. ”

Jiehoon Kim, MD, South Korea